





An easy-to-use medical image viewer for DICOM and certain non-DICOM medical images!

P.O. Box 8539 Greenville, SC 29604-8539 Tel: 864.979.1767 [www.vidistar.com](http://www.vidistar.com)


## Presenting ... VidiStar's Standalone Viewer™

- Applications include Emergency Medicine, Cardiovascular, General Ultrasound, Radiology, and other DICOM applications
- Microsoft XP and Vista compatible (Macintosh Compatibility Coming Soon!)
- Brightness/Contrast for still images & video clips (AVI, MPEG4, & DICOM)
- Open, review & compare multiple studies
- Drag & drop image(s) & video clip(s) for comparison




- DICOM and MPEG4 compatible
- Auto and manual calibration tools include measurement, pan/zoom, annotation, error correction, multiple image /video display, region of interest, free-hand surface area, magnify glass, and more

- Study List features include reports and study-related information
- Simply double click on the patient row to launch the viewer and view the study
- Synchronize and shuffle stress echo stages and views
- Burn CDs with VidiStar's View-Only Standalone Viewer for commercial burning applications



- Toggle view layout between 1:1, 2:1, 4:1, or 6:1
- Perform and calculate certain measurements (i.e Doppler)
- Compatible with existing systems such as SonoSite and other vendor DICOM modalities including Philips, GE, Biosound, ATL, and Siemens



- Easy to use graphical user-interface with point and click navigation
- Minimum computer requirements include at least 1 GB RAM and > 1 GHz processor

Contact VidiStar® today to purchase! 864.979.1767 or [sales@vidistar.com](mailto:sales@vidistar.com)